

DISCOVER TUBECLEAR

The Enteral Access Clearing System by Actuated Medical



Restore and maintain feeding tube patency without a trip to the emergency department.

Care for your patients without discomfort or unnecessary costs due to clogged feeding tubes.

Approximately 40% of TubeClear System inquiries are from patients or loved ones who are looking for a solution to minimize feeding tube replacement and keep the patient in their own bed.

When a clog cannot be cleared, patients may need to be transported to the hospital for another unclogging attempt or feeding tube replacement. Replacement renews risks of a misplaced feeding tube or aspiration pneumonia, increasing the emotional and monetary costs of care.

When used proactively, the TubeClear System keeps feeding tubes clear by removing sludge buildup along the inside walls of the tube.⁴ No need to transport to the hospital for unclogging or replacement, and no compromised enteral therapy delivery (i.e., medication, nutrition, hydration). Extended care patients on enteral nutrition reported 48 hours without nutrition due to a clogged feeding tube, increasing the risk of malnutrition.

How TubeClear Supports Extended Care Patients and Their Caretakers



Reduce need for hospital visits and transportation



Provide uninterrupted enteral therapy delivery



Offers comfortable bedside care to patients



Preserves facility HCAHPS reimbursements

Make every calorie count!



A **Clear** choice for extended care.

The TubeClear System is comprised of a single use Clearing Stem and a reusable Control Box. The Clearing Stem contains a flexible Wire surrounded by a Sheath. The Clearing Stem connects to the Control Box at the Diaphragm via the Magnet. The operator moves the Clearing Stem backward and forward using the motor within the Control Box, allowing the caregiver to adjust appropriately to different types of clogs.

TubeClear is indicated for use with most types of feeding and decompression tubes 6-18 French and 8-55 inches in length. Call us today to see how TubeClear can fit into your facility's goal of meaningful, cost-effective care.

**One clogged tube
can cost almost
\$12,000.**

How much can TubeClear save you?



No need to travel to the hospital for unclogging or tube replacement.

Potential Savings: \$460+



No missed medication or dehydration due to clogs.

Potential Savings: \$1000+



No risk of misplaced tubes during reinsertion procedures or associated infection.

Potential Savings: \$2000+



By improving patient care, your facility retains HCAHPS reimbursements while saving time and resources.

**Calculate your
potential savings at
[TubeClear.com/ROI](https://www.tubeclear.com/ROI)**

**Learn more at
[TubeClear.com](https://www.tubeclear.com) or call
814.355.003 ext 117
for more information.**

***This is where Actuated Medical can put the QR code to view sources, the TubeClear trademarked logo, additional logos/certifications.